

PROVOX[®] Dilator

**Rx
ONLY**

Provox Dilator REF 7211

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Contents

EN - ENGLISH	5
DE - DEUTSCH	8
NL - NEDERLANDS.....	11
FR - FRANÇAIS.....	14
IT - ITALIANO.....	17
ES - ESPAÑOL.....	20
PT - PORTUGUÊS	23
SV - SVENSKA.....	26
DA - DANSK.....	29
NO - NORSK.....	32
FI - SUOMI	35
IS - ÍSLENSKA	38
ET - EESTI.....	41
LT - LIETUVIŲ KALBA.....	44
CS - ČESKY.....	47
HU - MAGYAR.....	50
SK - SLOVENČINA.....	53
SL - SLOVENŠČINA.....	56
PL - POLSKI	59
RO - ROMÂNĂ.....	62
HR - HRVATSKI	65
SR - SRPSKI.....	68
EL - ΕΛΛΗΝΙΚΑ	71
TR - TÜRKÇE.....	75
HY - ՀԱՅԵՐԵՆ	78
RU - РУССКИЙ	82
MS - BAHASA MELAYU.....	86
JA - 日本語版.....	89
KO - 한국어.....	92
ZHTW - 繁體中文.....	95
HE - עברית.....	100
AR - العربية.....	103

Intended use

The Provox Dilator is intended for upsizing smaller tracheo-esophageal (TE) punctures (i.e. from 16 Fr diameter) to allow fitting of Provox voice prostheses, or upsizing a shrunken puncture to an adequate diameter i.e. after loss of a voice prosthesis.

The dilator may also be used for temporary blockage of a TE puncture or to temporarily prevent such from shrinkage.

It is only intended to be used by physicians or speech pathologists/therapists trained in the care and rehabilitation of laryngectomized patients.

Contraindications

- The device is not intended to be used for puncture dilation at the time of surgical creation of the puncture.
- Do not use the device in case of a small tracheostoma where the dilator could obstruct breathing.

Device description

Provox Dilator is a stepwise tapered, about 140 mm (5.5 inch) long solid curved rod made of medical grade silicone. The diameter is 15 Fr at the tip and increases to 24 Fr. At the end of each diameter step, i.e. 18, 20 and 22 Fr respectively, a small retaining collar is made to prevent the dilator from gliding back to the thinner section. The dilator also has a retainer strap with medallion intended to reduce the risk of accidental aspiration.

WARNINGS

- To avoid cross contamination, use appropriate hygienic procedures and make sure that each used device is cleaned, disinfected, dried and steam sterilized.
- Do not use the device if it shows any sign of structural damage.
- The material in the Provox Dilator withstands 50 cleaning and steam sterilization procedures at 137°C/278°F for 3 minutes or at 134°C/273°F for 18 minutes. Always ensure mechanical integrity before each use.
- Allergic reactions against the lubricating (anaesthetizing) gel may occur, read the instructions for use accompanying the lubricant. Always ask the patient about allergic history.

PRECAUTIONS

- Do not use any kind of oil or ointment for product care or during the procedure. Oil based lubricants (e.g. Vaseline) will weaken, damage or destroy the product and shall not be used at any time.
- Patients on anticoagulant therapy should be carefully evaluated for the risk of hemorrhage prior to dilation of the puncture.
- Make sure that the position chosen on the dilator is adequate for the puncture. Over-dilatation may result in leakage around the voice prosthesis.
- Never use excessive force to insert or remove the dilator from the puncture, this may cause bleeding or a false route.
- Do not use a contaminated dilator. Contaminated dilators may cause infection and inflammation.
- It is advisable to lubricate the dilator. Only use water-soluble lubricants in an appropriate amount.
- Do not use the dilator for a too long period of time as this may damage the esophagus wall.

Adverse events

Bleeding from the puncture during insertion may occur. Risk of bleeding can be reduced by inserting the dilator gently and slowly. Allow sufficient time for each section to dilate the puncture and do not force the dilator through the puncture to move to the next section.

Instructions for use

Checkout procedures

Always handle the dilator with clean gloves and avoid touching parts of it that will enter the puncture.

Check the integrity of the safety strap (3) by holding the device handle (1) and slightly pulling the medallion (2). Discard the device if it shows any sign of damage.

Prior to use, it is advisable to lubricate the dilator with some water soluble (anaesthetizing) gel.

Operating instruction

Remove the voice prosthesis from the TE puncture. Probe the puncture with the tapered tip of the dilator, its curvatures directed downwards, and gently insert the device until a certain puncture resistance can be felt.

Hold the dilator in this position for the desired amount of time, and then continue gently inserting it in the same way, section by section, until the desired puncture diameter has been reached. The waiting time needed for each section can vary from case to case. You shall carefully tape the medallion to the skin during the waiting time.

Once the dilator has been in the final position for the desired amount of time, gently pull the rod - not the medallion - for removal and insert the Provox voice prosthesis using its standard procedure.

Cleaning and sterilization

Cleaning, disinfection and sterilization are described in “Instructions for cleaning and sterilization”. If this information is missing, please contact your supplier or visit our website on www.atosmedical.com.

The material in the Provox Dilator withstands 50 cleaning and steam sterilization procedures at 137°C/278°F for 3 minutes or at 134°C/273°F for 18 minutes.

How supplied

The device is manufactured in a cleanroom and supplied non-sterile.

Device lifetime

Caution: The dilator shall be replaced every 2 years, or after 50 sterilization cycles, whichever comes first. Check mechanical integrity of the device prior to each use according to the Checkout procedures.

Disposal

Always follow medical practice and national requirements regarding biohazards when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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